CooperVision®

Please find the Instructions for Use (IFU) for the following products registered with Medical Device Authority in Malaysia.
Sila rujuk kepada Arahan Penggunaan bagi produk berikut yang berdaftar dengan Pihak Berkuasa Peranti Perubatan di Malaysia.

NO. Bil.	PRODUCT NAME Nama Produk	MDA REGISTRATION NUMBER Nombor Pendaftaran MDA
1	ASCEND 1 DAY AIR (SOMOFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB81158769718
2	ASCEND PREMIER (COMFILCON A) MONTHLY CONTACT LENSES	GB68048769918
3	BIOFINITY (COMFILCON A) MONTHLY CONTACT LENSES	GB13979771318
4	BIOFINITY MULTIFOCAL (COMFILCON A) MONTHLY CONTACT LENSES	GB25597771618
5	BIOFINITY TORIC (COMFILCON A) MONTHLY CONTACT LENSES	GB37084771518
6	BIOMEDICS 1 DAY EXTRA (OCUFILCON D) DAILY DISPOSABLE CONTACT LENSES	GB24557640718
7	BIOMEDICS 1 DAY EXTRA TORIC (OCUFILCON D) DAILY DISPOSABLE CONTACT LENSES	GB13426771718
8	BIOMEDICS 55 EVOLUTION (OCUFILCON D) MONTHLY CONTACT LENSES	GB39364639818
9	BIOMEDICS NOW (OCUFILCON D) MONTHLY CONTACT LENSES	GB33494769618
10	BIOMEDICS TORIC (OCUFILCON D TORIC) MONTHLY CONTACT LENSES	GB77143771418
11	CLARITI 1 DAY (SOMOFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB98072772318
12	CLARITI 1 DAY MULTIFOCAL (SOMOFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB99714639618
13	CLARITI 1 DAY TORIC (SOMOFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB97189772418
14	CLARITI ELITE (SOMOFILCON A) MONTHLY CONTACT LENSES	GB42743483117
15	EYECONT PRO COMFORT (OMAFILCON B) MONTHLY CONTACT LENSES	GB81362771818
16	EYECONT PRO COMFORT DAILY (OCUFILCON D) DAILY DISPOSABLE CONTACT LENSES	GB19110819818
17	LIVE 1 DAY (SOMOFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB8742020-42389
18	MISIGHT 1 DAY (OMAFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB10878019-37149
19	MYDAY (STENFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB93422640618
20	OCUFILCON D 55% 1 DAY CONTACT LENSES	GB65158639518

NO. Bil.	PRODUCT NAME Nama Produk	MDA REGISTRATION NUMBER Nombor Pendaftaran MDA
21	PRO-FIT PERFORMANCE + (FANFILCON A) MONTHLY CONTACT LENSES	GB181471020518
22	PRO-FIT PERFORMANCE + TORIC (FANFILCON A TORIC) MONTHLY CONTACT LENSES	GB83150771218
23	PROVISION 55 ASPHERIC (OCUFILCON D) MONTHLY CONTACT LENSES	GB53428771918
24	PROVISION 55UV (OCUFILCON D) MONTHLY CONTACT LENSES	GB24429640018
25	SIMVUE BIO MOIST 1 DAY (OCUFILCON D) DAILY DISPOSABLE CONTACT LENSES	GB31003939518
26	SIMVUE BIO MOIST 55 (OCUFILCON D) MONTHLY CONTACT LENSES	GB97310771118
27	SMARTVIEW AQUA (COMFILCON A) MONTHLY CONTACT LENSES	GB16642770518
28	SMARTVIEW AQUA ONE DAY (SOMOFILCON A) DAILY DISPOSABLE CONTACT LENSES	GB43196770718
29	SOFTVIEW 1 DAY EXTRA (OCUFILCON D) DAILY DISPOSABLE CONTACT LENSES	GB44637640118
30	SOFTVIEW 55 UV (OCUFILCON D) MONTHLY CONTACT LENSES	GB23838772018
31	SOMOFILCON A DAILY DISPOSABLE CONTACT LENSES	GB55065483217
32	ULTRAFLEX II 55UV (OCUFILCON D) MONTHLY CONTACT LENSES	GB13278772118
33	ULTRAFLEX PLUS (SOMOFILCON A) MONTHLY CONTACT LENSES	GB620601123918

SOFT CONTACT LENSES

Instructions for Use

A Guide to Trouble-Free Soft Contact Lens Wear

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Local Authorised Representative:

Cooper Vision Contact Lens Malaysia Sdn. Bhd.
No. 50-1-11 & 12, Level 1, Wisma UOA Damansara,
No. 50, Jalan Dungun, Damansara Heights,
50490 Kuala Lumpur, Malaysia.

www.coopervision.com.my

Always wash and dry your hands before handling your lenses.

Do not use if packaging is not sealed.

Contact lenses are not to be shared.

LENS INSERTION

1. Before insertion, ensure the lens is not inside out. Place on your index finger to verify the shape (Fig. 1 & 2).
2. Using the middle finger, hold your lower lid down. With the index finger of the other hand, raise upper lid. Place the lens onto the eye (Fig. 3). Repeat for the other eye.

LENS REMOVAL

Separate the lids with the middle finger of each hand. Look up and with the index finger of the lower hand; slide the lens down onto the white of the eye. Pinch the lens and lift out with the thumb and index finger (Fig. 4). Repeat for the other eye.

DISINFECTION/SOAKING

For lenses worn more than once clean and disinfect them following strictly the instruction on your solution bottle or as directed by your eyecare practitioner.

Daily disposable lenses are intended to be worn only once and then discarded. If your lenses are intended for daily disposable wear, **DO NOT** wear them more than once. Daily disposable lenses are not designed to be cleaned or disinfected and worn more than one time. Daily disposable lenses worn more than one time may become damaged and could put you at great risk of contact lens care and handling related eye infections or other problems.

IF THE LENSES DRY OUT:

Re-wet with disinfecting, soaking or rinsing solution then follow the Disinfecting/Soaking step above.

PRECAUTIONS:

1. For all types of lenses, your eyecare practitioner will recommend a lens wearing, cleaning and replacement regime. Do not vary these. If in doubt seek the advice of your eyecare practitioner.
2. If wearing lenses on an extended wear basis, remove and dispose of your lenses the evening before you are due to replace them, inserting new lenses the following day.
3. For lenses worn more than once, always clean and disinfect your lenses as instructed, after lens removal.
4. If your eyes become red or irritated or you experience any pain or unexpected change in vision, remove your lenses immediately and consult your eyecare practitioner.

5. Avoid contaminating lenses with make-up, aerosol sprays or wearing them near noxious fumes.
6. Replace damaged, contaminated or scratched lenses immediately.
7. Avoid long nails, which harbour bacteria and could damage the eye or lens. A scratched eye is vulnerable to infection.
8. Keep your lenses out of the reach of children.
9. Only use your lenses for water sports if watertight goggles are worn.

WARNINGS:

- 1) **DO NOT** sleep in your lenses unless your eyecare practitioner has advised it is safe to do so.
- 2) **DO NOT** wear your lenses beyond the period recommended by your eyecare practitioner.
- 3) **DO NOT** use household products (e.g. disinfectants) on your lenses.
- 4) **DO NOT** wear your lenses when using eye drops (unless advised to do so by a medical practitioner).
- 5) **DO NOT** use saliva to wet your lenses.
- 6) **DO NOT** rinse your lenses or lens case with tap water.

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